



Bulb pipettes, made of plastic, 1 mark

PP, high clarity. Calibrated to deliver (TD, Ex). Resistant to breakage. Exposure to temperatures up to 60 °C will not cause permanent exceeding of error limits! Cleaning temperature below 60 °C is recommended to preserve marks and inscriptions. 1 and 2 ml: pack of 12; 5 to 50 ml: pack of 6.

Capacity ml	Pattern	Error limit ± ml	Length mm ± 10 mm	Cat. No.
1	straight	0.02	300	300 02
2	straight	0.02	300	300 03
5	bulb	0.03	300	300 07
10	bulb	0.04	440	300 12
25	bulb	0.06	450	300 15
50	bulb	0.1	460	300 18